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Aerosol a micro pistone
Micro Piston Nebuliser
Aérosol à micro-piston
Mikrokolben-aerosol
Aerosol de micropistón
Aerossol de micro pistão
Inhalator mikrokompresorowy
Inhalator met micro zuiger
البخاخة الكباسية الصغيرة

REF 00010230000100

FIGURA A/ FIGURE A/ ABBILDUNG A/ RYCINA A/ FIGUUR A / A - الشكل



FIGURA / FIGURE / ABBILDUNG / RYCINA / FIGUUR / الشكل

ISTRUZIONI D'USO

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EN USER MANUAL

Dear customer, thank you for choosing this Chicco Micro Piston Nebuliser, a nebuliser therapy unit designed and manufactured according to the most up-to-date technologies. The nebuliser is a class IIa medical device according to European Regulation for medical devices EU Reg 2017/745. Before using the unit, we highly recommend that you read the shortlist of warnings contained in this manual, to make sure you have correctly understood how to use it.

IMPORTANT:

Read the instructions carefully before use.

COMPONENTS (PICTURE A):

1. Compressor Housing - 2. Air Outlet - 3. Filter Compartment - 4. Usb C Port - 5. On/Off Button- 6. Nebuliser cup with linear connector - 7. Mouthpiece - 8. Nasal Fork - 9. Air Tube - 10. Adult Mask - 11. Paediatric Mask - 12. 5 Replacement Filters - 13. Type-C Usb Cable - 14. Elbow Connector.

INTENDED PURPOSE:

The device is an aerosol therapy system suitable for domestic use.

The device is designed for the production of compressed air to operate a small volume nebuliser for the production of medical aerosols for lung respiratory disorders.

PATIENT POPULATION:

The device is intended for use with children from 2 years old, adolescent and adult patients.

INTENDED USERS:

The user is an adult person without specific knowledge or professional ability. The patient is the intended operator except in the case of children and patients that requires special assistance. If children or people with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device.

INDICATIONS:

Lung acute or chronic diseases of respiratory disorders organs, or inflammation of the upper respiratory system.

CONTRAINDICATIONS:

The device is not indicated to use with quick-relief medications during life-threatening asthma attacks. No contraindications exist to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

WARNINGS:

Please report any serious incident that has occurred in relation to the device injury or adverse event to the competent local authority and to the Manufacturer or to the European Authorized Representative (EC REP) Vigilance contact points: https://health.ec.europa.eu/medical-devicessector/new-regulations/contacts_en

- Use the device only as described in this manual and therefore as an aerosol therapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- Keep the user manual for future reference.
- Do not operate the unit in the presence of any anesthetic mixture inflammable with oxygen or nitrogen protoxide.
- The correct functioning of the equipment can be affected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it into a different power socket.
- In case of failure and malfunction, read the "Malfunctions and Actions to take" in the user manual. Do not handle or open the compressor housing.
- For repair operations, address only a technical service center authorized by the manufacturer and require the use of original spare parts. The non-observation of the above-mentioned indications can compromise the device's safety.
- Comply with the safety regulations concerning electrical devices and in particular:
 - use only original applied parts;
 - never submerge the unit in water;
 - never wet the device, it is not protected against water penetration;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - place the unit on a stable and horizontal surface during its operation;
 - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
 - the power plug is the separation element from the grid power; keep the plug accessible when the device is in use.
- Before plugging into the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated or has fully charged.
- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause

- damage to people, animals or things, for which the manufacturer cannot be held responsible.
- The USB cable cannot be replaced by the user. In case of damage, address to a technical service center authorized by the manufacturer for its replacement.
- The power supply cord should always be fully unwound in order to prevent dangerous overheating.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply.
- Some parts of the unit are so small that they may be swallowed by children; keep the equipment out from children's reach.
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
 - use this device only with medicines prescribed by your doctor.
 - Only administer the treatment using the applied part recommended by your doctor depending on the pathology.
 - Use the nasal fork only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose but only bringing them as close as possible.
- Check the medicine package leaflet for possible contraindications for use with common aerosol therapy systems.
- To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
- Do not position the equipment so that it is difficult to operate the disconnection of the device.
- Nebuliser and applied parts are single patient reuse. Device is multi-patient reuse.
- Nebulizing system is not suitable for use in anesthetic breathing system or a ventilator breathing system. The patient cannot undergo MRI scan while using this device.
- Device is designed to nebulize solution and suspension liquids
- Do not use the device if you spot damage or notice something unusual.
- This device consists of delicate components and must be handled with care. Observe the storage and operating conditions in the chapter "Technical Data"
- The use of this device does not replace a visit to the doctor. If there is no improvement in your health after treatment, consult your doctor again.

HOW TO USE THE DEVICE

The device must be checked before each use in order to identify any malfunctions and/or damage caused during transportation and/or storage. During inhalation, sit upright and relaxed to avoid compressing the airways and impairing treatment effectiveness. The applied parts must be used only with a single patient, it is not recommended to use them with several patients.

- After unpacking the device, check it for visible damage or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check if the applied parts are intact.
- Before using the device, proceed with the cleaning operations as described in the chapter "CLEANING AND MAINTENANCE". Applied parts can resist dishwasher and steam treatment.
- Open the nebuliser by turning the top counterclockwise (**PICTURE B**).
- Make sure that the drug delivery cone is inside the nebulizer ampoule and that the linear connector is assembled onto the ampoule. (**PICTURE C**).
- Pour into the nebuliser the amount of medicine indicated by your doctor without exceeding the maximum level of the nebulizer cup (**PICTURE D**).
- Close the nebuliser by turning clockwise the two parts making sure that they are properly sealed (**PICTURE E**).
- Connect one end of the air tube to the dedicated outlet at the bottom of the nebuliser (**PICTURE F**) and the other end to the air outlet (**PICTURE G**) on the device.
- Use elbow connect To connect the mask to the ampoule using the linear connector. The nasal fork and mouthpiece can be connected directly to the ampoule via the linear connector. (**PICTURE H**)
- To start treatment, connect the device to the chosen power source via the USB cable, paying attention to compliance with the power specifications, as indicated in the TECHNICAL SPECIFICATIONS paragraph. (**PICTURE I**)
- Set the I/O switch into the "I" position. (**PICTURE J**)

CAUTION:

**The device is designed for intermittent use (30 min. on - 30 min. off)
Do not use the product continuously for 30 minutes or longer. Using the device 30 minutes before it have elapsed may result in motor overheating and cause the safety cut-off devices to activate. To prevent overheating and damage to the compressor, switch off the appliance for at least 30 minutes after every 30 minutes of uninterrupted use.**

- Inhale the aerosol solution using the prescribed applied part.
- Once you have finished the treatment, switch the appliance off by setting the I/O into the "O" position and disconnect the plug from the wall socket.
- Clean the nebuliser and its applied parts following the steps in the «CLEANING AND MAINTENANCE» section.
- The device requires no calibration.
- Tampering with the device is absolutely forbidden.
- No modification to the device is permitted.

CLEANING AND DISINFECTION

DEVICE(S): AEROSOLTHERAPY APPLIED PARTS	
WARNINGS	Before the first use and after each treatment. Follow carefully the cleaning and disinfection instructions of the applied parts as they are very important for the device performances and the therapy success. Use original parts only. Do not clean or disinfect the air tube. The applied parts can not be cleaned and disinfected by automated method.
Limitations on processing	The nebulizer must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebulizer nozzle is obstructed by dry medicine, dust, etc. We recommend to replace the nebulizer after a period between 6 months and 1 year, depending on the usage. The maximum times of cleaning and disinfection of the nebulizer are 360 times. The maximum times of cleaning and disinfection of the mask, the nosepiece and the mouthpiece are 360 times.

INSTRUCTIONS	
Preparation before cleaning	Detach the air tube from the nebulizer. Detach the mouthpiece or nosepiece or mask from the nebulizer. Make sure that all residual volume is removed from the nebulizer. Disassemble the nebulizer by turning counterclockwise the top and remove the medicine conduction cone. Note: To avoid microbial growth and residual drug drying, clean and disinfect the applied parts immediately after each use.
Cleaning	Briefly rinse all parts in running tap water for at least 10 seconds beforehand. Mix little dish detergent and warm tap water (i.e. FAIRY brand dish washing liquid for hand washing in a ratio of 2ml: 1L) in a clean container. Immerse the components of disassembled nebulizer, the mask, the mouthpiece and the nosepiece into the mixed water for about 5 min. And then scrub all the surface of all the components with clean and small brush for at least 8 times. Afterwards rinse all parts thoroughly in running tap water at least 30 seconds, to fully remove any possible dish washing liquid residue.
Disinfection	After cleaning, disinfect all of the disassembled parts (only parts that have been cleaned can be disinfected effectively). Boil the disassembled nebulizer, the mask, the mouthpiece and the nosepiece for 5 minutes in boiling tap water.
Drying	Reassemble the nebulizer components and connect it to the air-outlet, switch the device on and let it work for 10-15 minutes. Let all parts dry completely before re-assembling and reuse to avoid risk of germ grow. Care needs to be taken not to contaminate the parts after they have been cleaned and disinfected. Contamination can be avoided by good hand washing and not touching the inside sections of the device when laying them out to dry or when reassembling.
Inspection	Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts.
Packaging	Pack dry parts in a clean and sealed container when not in use. Do NOT pack wet or damp parts.
Storage	Storage conditions refer to "TECHNICAL SPECIFICATIONS". Note: re-clean and re-disinfect the parts if they are stored more than one day.
Transportation	After cleaning and disinfection, always transport the parts with clean and sealed container. Contamination can be avoided by good hand washing and not touching the inside sections of the parts when take out and re-assemble the parts for use.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. The processor retains the responsibility to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

FILTER REPLACEMENT

Periodically check the condition of the filter. The filter was inserted to protect the compressor. Correct filter maintenance prolongs the life of the nebuliser. The frequency with which filter should be replaced depends on the conditions in which the device is used. In normal usage conditions, the air filter must be replaced approximately after 100 working hours or once a year. We recommend to periodically check the air filter (10-12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Remove the filter and replace it with a new one. Do not attempt to clean the filter to reuse it. The air filter shall not be serviced or maintained while in use with a patient.

 **Only use original filters. Do not use the device without a filter.**

MALFUNCTIONS and ACTIONS TO TAKE

PROBLEM	POSSIBLE SOLUTION
The device fails to switch on	• Make sure the plug is properly inserted in the power outlet. • Ensure the ON/OFF switch is in the position «I». • Make sure that the device has been operating within operating limits indicated in this manual (30min ON/30min OFF). • Make sure that the power supply meets the requirements specified in Technical Specification (5V >2A output). If the problem persists it is suggested to change power supply.
The device does not nebulize or nebulizes weakly.	• Make sure the air hose fittings are properly inserted in the compressor and the nebuliser connections. • Ensure the air tube is not squashed, bent, dirty or blocked. If necessary, replace it with a new one • Make sure the nebuliser is not empty or that it has been filled with the correct amount of medication (MAX. 6 ml). • Make sure that the nebuliser nozzle is not clogged.
The device is noisy	• Make sure that the filter is mounted properly
Low nebulisation	• The nebuliser cup can be clogged. Remove and detach all the nebuliser cup components and clean them as illustrated in this instruction manual. If the issue persists after cleaning, replace the nebuliser cup.
No nebulisation	• Nebuliser cup components can be not assembled correctly. Disassemble and reassemble the nebuliser cup; if the issue persists, replace the nebuliser cup. • Check that the device is connected to the chosen power source via the USB cable and that the on/ off switch is on. • If the issue persists, contact customer service.
Slow nebulisation	• Drug may be too thick. Dilute the drug as instructed by the manufacturer.

In the event that the equipment does not resume correctly, please **contact CUSTOMER SERVICE**.

MAINTENANCE AND REPAIRS

In case of failure, contact **Customer Service**. **Never open or try to repair the device**. The device has no user-serviceable parts within and does not need internal maintenance or lubrication.

TECHNICAL SPECIFICATIONS

Item	Specification	
Max Pressure	70-120	kPa
Operating Pressure (with calibrated nozzle, Φ 0,50)	35-60	kPa
Protection type IP	IP 21	--
Sound noise level (Distance: 1m, Front direction) max level	≤ 48	dBA
Equipment conditions of use	30 ON / 30 OFF	min
Dimensions	106x58x45	mm
Weight	0,17 $\pm$ 0,05Kg	kg
Service life	400	Ore
Power input	5VDC > 2A	--
Current ABSORPTION	≤ 2000	ml/min
Voltage power source	5VDC	V
Neb Rate (average value)	0,3	ml/min
Neb Rate (acceptable range)	0.2-0.5	ml/min
Nebuliser kit AEROSOL PERFORMANCES ACCORDING TO EN ISO 27427:2019 based on adult ventilatory pattern with NaF		--

Aerosol output	0,51	ml
Aerosol output rate	0,16	ml/min
Percentage of fill volume emitted per min	8	%
Residual Volume	1,47	ml
Particle size (MMAD)	2,9	µm
GSD (geometric standard deviation)	2,8	- -
RF (respirable fraction <5µm)	71,1	%
Large particle range (>5µm)	29	%
Middle particle range (2 to 5µm)	26	%
Small particle range (<2µm)	45	%
Min & Max filling volume:	2 - 6	ml
Operating Temperature	10 ~ 40	°C
Operating Humidity	10 ~ 95	% RH
Operating Pressure	700 ~ 1060	hPa
Storage Temperature	-20 ~ 60	°C
Storage Humidity	10 ~ 95	% RH
Storage Pressure	700 ~ 1060	hPa
Protection against electrical shock	2	/

- Class II device referred to protection against electric shock.
- Nebuliser, mouthpiece, masks and nasal fork are type BF applied parts.
- Important for the safety of the device and its user: use only power supplies as specified (5V >2A output with USB-C connector and ensure their certification in accordance with EN 60601-1 and EN 60601-1-2)

Note: Technical specifications subject to change without notice.

AEROSOL CHARACTERISTICS IN ACCORDANCE WITH REGULATION EN ISO27427

Standards applied:

Electric safety standards CEI EN 60601-1 Electromagnetic compatibility according to CEI EN 60601-1-2 The device is a Class IIa medical device. **Device Complies with European Regulation on Medical Devices EU MDR 2017/745.**

DISPOSAL:

ITEM	COMPONENT	DISPOSAL DESCRIPTION	PHOTO
1	Device	The component mainly includes plastic and electronic components. All comply with RoHS and REACH standards and are suitable for safe disposal. This product is subject to European Directive 2012/19/ EU on electrical and electronic equipment waste and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products.	
2	Applied parts (tube, nebuliser cup, elbow connector, mouthpiece, nasal fork, mask)	The components are made of plastic. All comply with RoHS and REACH. Before disposal, all the parts should be cleaned according to manual instructions and then disinfected by boiling for 5 minutes.	

WEEE- Waste Electrical and Electronic Equipment

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions. The warranty shall not therefore apply in case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national laws in the country of purchase, where appropriate.

DEVICE SYMBOLS - MANUFACTURER'S DECLARATION AND GUIDELINES - ELECTROMAGNETIC EMISSIONS: SEE APPENDIX.

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DEVICE SYMBOLS

	Reference Number		Distributor
	ON /OFF		Authorized representative
	Alternating current		Manufacturer
	Class II		Date of manufacture (year, month and day)
	Type BF		Serial Number (LOT/serial number)
	Device Complies with European Regulation on Medical Devices EU MDR 2017/745		Caution
	Medical device		Humidity limitation
	Follow instruction for use		Temperature limitation
	Importer		Ambient pressure limitation
	Unique device identifier		Type Number
	Batch code		Single patient multiple use (for applied parts only)
	Protection against solid foreign objects and against harmful effects due to the ingress of water		

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic emissions			
The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 2\text{kV} \pm 4\text{kV} \pm 6\text{kV} \pm 8\text{kV} \pm 15\text{kV}$	$\pm 8\text{ kV}$ contact $\pm 2\text{kV} \pm 4\text{kV} \pm 8\text{kV} \pm 15\text{kV}$ air	The floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	$\pm 1\text{kV}$, $\pm 2\text{kV}$, 100 kHz repetition frequency	Power supply lines: $\pm 2\text{ kV}$	The main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0.5\text{kV}$, $\pm 1\text{kV}$ (Line to line) $\pm 0.5\text{kV}$, $\pm 1\text{kV}$, $\pm 2\text{kV}$ (Line to Ground)	Line to line: $\pm 0.5\text{kV}$, $\pm 1\text{kV}$ Line to Ground: $\pm 0.5\text{kV}$, $\pm 1\text{kV}$, $\pm 2\text{kV}$	The main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T for 0.5 cycle 0% U_T for 1 cycle 70% U_T for 25/30 cycles 0% U_T for 250/300 cycles	0% U_T for 0.5 cycle 0% U_T for 1 cycle 70% U_T for 25/30 cycles 0% U_T for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) Magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic emissions			
The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Smart including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. Recommended separation distances $d = 1,2 \sqrt{P}$ $d = 2 \sqrt{P}$

Radiated RF IEC61000-4-3	10V/m	10V/m	80 MHz a 800MHz $d=1.2 \sqrt{P}$ Da 800 MHz a 2.7 GHz $d=2.3 \sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the device

The Device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d=1.2\sqrt{P}$	150 kHz to 80 MHz (in ISM and amateur radio bands) $d=2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.7GHz $d=2.3\sqrt{P}$
0,01	0.12	0.2	0.12	0.23
0,1	0.38	0.632	0.38	0.73
1	1.2	2	1.2	2.3
10	3.8	6.32	3.8	7.3
100	12	20	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

	Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{a)}	Modulation ^{b)} (W)	Distance (m)	Immunity Test Level (V/m)
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
	450	380 – 390	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
	870						
	930						
	1 720	1 700 – 1 990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
	1 845						
	1 970						
	2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
	5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
	5 500						
	5 785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50% duty cycle square wave signal.

^{c)} As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The Manufacturer should consider reducing the minimum separation distance, based on Risk Management, and using higher Immunity Test Levels that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher Immunity Test Levels shall be calculated using the following equation:

$$E = 6/d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the Immunity Test Level in V/m

NOTE

[illegible]

NOTE

This image shows a full page of white paper with horizontal blue dashed lines, typical of primary school writing paper. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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REF 00010230000100

Type: GCE863

Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road,
European Industrial Zokeneo quellne, Xiaolan Town, 528415 Zhongshan City,
Guangdong Province, PEOPLE'S REPUBLIC OF CHINA,
Phone: +86 760 22589901

Donawa Lifescience
Piazza Albania, 10 ; 00153 Roma Italia
<https://www.donawa.com/wli/main/contatti.index>

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Artsana S.p.A