Aerosol a Pistone
SUPER SOFT

Super Soft Nebulizer

REF 00009067000000

REF 00009813000000

Fig. 1


Fig. 2


Fig. 3




fig. 3C

ISTRUZIONI D'USO

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Numero verde consumatori Italia 800 188 898

EN INSTRUCTIONS

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,

Thank you for choosing this Chicco aerosol piston, an aerosol therapy unit designed and manufactured according to the most up-to-date technologies.

The aerosol is a class IIa medical device which converts medicines from liquid or suspension form to aerosol form, so they can be administered directly to the airways. The nebuliser cup supplied with the product allows the medicine to correctly penetrate the respiratory tract, thanks to size of the aerosol particles.

Additionally there is a Nasal irrigation device **REF** 00009066000000, ideal for washing nasal cavities. The nasal irrigation device is compatible with Chicco nebulizers Super Soft **REF** 00009067000000 and **REF** 00009813000000 and may be purchased separately.

Before using the unit, we highly recommend that you read the shortlist of warnings contained in this manual, to make sure you have correctly understood how to use it.

ACCESSORIES

A - Adult mask

B - Children's mask

C - Mask connector

C1 - Mouthpiece

D - Nose fork

E - Connection tube

F - Nebuliser cup

A, B, C, C1, D, E, F are available as spare parts

in kit **REF** 00009069000000

DEVICE COMPONENTS

G - On/Off switch

H - Compressed air outlet

I - Filter

L - Filter cover cap

M - Nebuliser cup support

N - Openings for cooling

O - Power cord

P - Item holder compartment lid

Q - Item holder compartment

R - Lid Baby **REF** 00009067000000 / **REF** 00009813000000

GENERAL WARNINGS

- Before use, carefully read the information in this manual, and keep for future reference.
- Unit designed for aerosol therapy. Follow the instructions of a medical practitioner as to the type of medicine to use, the dose, frequency and duration of inhalations. Any use other than the one for which the unit is intended is considered improper and therefore dangerous; Artsana cannot be held liable for any damage caused by incorrect or unreasonable use or if the device is used in electrical systems that do not conform to safety regulations in effect.

- A power outage, a sudden failure or any other adverse condition may cause the unit to stop working, therefore it is recommended that users have an alternative unit or medicine (agreed upon with a physician) that can be used in the event this occurs.
- After removing the device from its packaging, check that it is intact, without any visible damage that may have occurred during transport. If in doubt, do not use the unit and contact Customer Services.
- Unroll the power cord completely before connecting it to the mains.
- The incomplete unrolling of the power cord may cause overheating.
- Packaging (bags, box, etc.) may be dangerous and should be kept out of reach of children. Keep the unit out of reach of children, pets and persons with cognitive limitations.
- Before connecting the device, check that your mains voltage corresponds to that shown on the rating label. The rating label is located on the bottom of the unit.

In case of incompatibility between the unit plug and the electrical socket, use certified adapters if the existing legislation in the individual country allows for it or have an Artsana S.p.A. authorised technician or Customer Services replace the plug.

To prevent overheating and damage to the compressor, switch off the appliance for at least 40 minutes after every 20 minutes of uninterrupted use.

• WARNING! Using the device before 40 minutes have elapsed may result in motor overheating and cause the safety cut-off devices to activate.

• To ensure correct operation, make sure that the air filter is dry. The filter cover may be removed.

• WARNING! keep the filter cover and the filter out of the reach of children as these are small parts and could represent a swallowing and suffocating hazard.

• The Baby (R) lid is not a toy: keep out of children's reach.

- Do not use the device in the presence of nitrous oxide, oxygen or flammable air/anaesthetic mixtures.
- Turn off the device and disconnect from the mains after use and before adding more medicine. Do not fill the nebuliser cup beyond the maximum level indicated on it and check that it is assembled with all its components.
- Keep the device and power cord away from hot surfaces.
- Do not use the device while taking a bath or shower, or when in a damp place, or close to bath-tubs, sinks, washbasins, or in any other situation in which there are liquids that may come into contact with the device.
- Never touch the appliance with wet or damp hands.
- Do not drop or put the device into water or other liquids. Should this occur, unplug it immediately, stop using the appliance and contact Artsana S.p.A. or Customer Services.
- Do not block the openings for cooling air.
- Do not use the device if you feel drowsy.
- If children or persons with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device. The device contains small parts that can be swallowed and the power cord poses risk of strangling.
- Do not use connectors or accessories not recommended by the manufacturer.
- Make sure the unit is placed on a level, stable surface when used in order to prevent spillage.
- During use, the surface must be free of any objects that may obstruct the correct flow of air.
- Do not block the holes of the filter cover located on the lower surface of the device with cloth

- or other items: this could cause the motor to overheat and reduce device performance
- Do not pour liquid on the device, especially on the filter cover: this could lead to the undesired discharge/nebulisation of he liquid.
 - Before cleaning or performing any maintenance work on the device, disconnect it from the mains by pulling the plug out of the electrical socket.
 - If you decide to no longer use the device, disconnect from the mains and cut off the power cord so the device cannot work. Dispose of the device and power cord in compliance with existing laws. We also recommend eliminating any parts that are potential hazards, particularly for children.
 - Refer to laws in force for information regarding the disposal of accessories subject to wear. To dispose of the device, refer to Directive 2012/19/EU.
 - Do not use the device if, after a fall, it shows signs of damage on any of its parts. If in doubt contact Customer Services.
 - In the case of the failure or malfunction of the appliance, switch it off by removing the plug from the mains and contact Customer Services only. The power cord and fuse specified in the "Technical Specifications" must only be replaced by Customer Services staff.
 - **WARNING! The power cable is the separating element from the electrical mains. When in use, to disconnect the appliance from the mains, switch it off and pull out the plug from the power outlet. Always connect the plug only to a power socket that is easily reachable to facilitate disconnection from the mains. Do not attempt to open, repair or tamper with the appliance.**

PREPARING AND USING THE APPLIANCE (fig. 3)

- Take the device and the accessories out of the box.
- Pour the medicine and/or saline solution in the nebuliser tank, according to the doses recommended by your physician, without exceeding the maximum level shown by the mark on the tank; close the nebuliser cup;
- Connect the nebuliser cup to the accessory needed for the therapy: Adult or children's mask (attach the mask connector to the mask you intend to use fig. 3C), fork, mouthpiece (applied parts). Ask your doctor for suggestions on which accessory to use if you are not sure;
- Connect the tube to the compressed air outlet and to the nebuliser cup;
- Connect the device to the mains;
- Turn on by pressing the on/off switch (G) and start therapy;
- Carry out treatment preferably seated and in a comfortable position;
- When the aerosol jet becomes intermittent, interrupt the treatment for a few seconds to allow the suspended drops of solution to fall on the nebuliser cup walls. Resume the treatment session, which will be finished when no more nebulised solution is released from the nebuliser cup;
- Once the treatment session is completed, switch the device off by unplugging it, remove the accessory used and clean the device and accessories according to the instructions in the Cleaning and Maintenance paragraph;
- The patient and operator are the same person. The steps that the patient/user can perform safely are only those specified in the paragraphs "Preparing and using the appliance," "Cleaning and maintenance of the compressor", "Cleaning and maintenance of the nebuliser cup and accessories." Warning! Do not service and/or repair the appliance while it is being used! Any other maintenance/service activity/operation must be carried out solely by Customer Services.

COMPRESSOR CLEANING AND MAINTENANCE

- **WARNING! Do not wet the device or put it into water or other liquids. For cleaning, use a clean, dry cloth only.**

CLEANING AND MAINTENANCE OF THE NEBULISER CUP AND ACCESSORIES

- Remove and detach all the nebuliser cup components and dip them in boiling water for about 5 minutes.
- **WARNING! Do not boil the masks.**
- Rinse the accessories (masks, connector, mouthpiece, fork and tube) under running water. Cold disinfect all parts that come in contact with the patient using denatured alcohol. Carefully dry all accessories and the device before reassembling. Store the device and the accessories in a cool, dry place, away from light and heat.
- Never use benzene, thinners or other inflammable chemical substances for cleaning.
- For increased hygienic safety, never use the same accessories on more than one patient. Purchase a dedicated kit for each user (Nebuliser Kit Super Soft **REF** 00009069000000).
- Clean and check that the condition of the accessories before and after use. Replace them if damaged.

CHECKING AND REPLACING FILTERS

Periodically check the condition of the filter. The filter was inserted to protect the compressor. Correct filter maintenance prolongs the life of the aerosol. The frequency with which filter should be replaced depends on the conditions in which the device is used. The filter is located on the back of the device.

TO REPLACE THE FILTER (fig. 2):

- Remove the filter cover cap turning it anti-clockwise;
- Remove the filter to be replaced from its housing (using a small tool if necessary);
- Carefully fit the new filter in its housing;
- Replace the filter cover cap and fasten tightly by turning it in a clockwise direction;
- **WARNING! Never leave the removed filter cover cap and/or filters unattended: These are small parts posing the risk of suffocation if swallowed by children.**

TYPE OF DEFECT	CAUSE	REMEDY
Scarce nebulisation	Clogged nebuliser cup	Remove and detach all the nebuliser cup components and clean them as illustrated in this instruction manual.
Scarce nebulisation	Clogged nebuliser cup	If the issue persists after cleaning, replace the nebuliser cup.
No nebulisation	Nebuliser cup components not assembled correctly	Disassemble and reassemble the nebuliser cup; if the issue persists, replace the nebuliser cup.
No nebulisation	The compressor does not work	Check that the power cord plug is disconnected from the mains and that the on/off switch is on. If the issue persists contact Customer Services.
Slow nebulisation	Drug too thick	Dilute the drug as instructed by the manufacturer.

TECHNICAL DATA

Nominal voltage: 230 V~

Service life in storage conditions: 5 years

Absorbed current: 0.6 A

Frequency: 50 Hz

Maximum flow: 14 ± 4 l/min

Fuse: F5 L 250V, 5x20mm

Max pressure: 2.55 ± 0.6 bar

Operational flow: 6 ± 2 l/min

Conditions for using the device: 20' ON - 40' OFF at 40 °C

Operational pressure: 0.76 ± 0.21 bar

MMAD: $< 3.7 \mu\text{m}$

Weight: 1.6 kg

Nebuliser cup capacity: 10 ml Max; 2ml Min

Breathable fraction: $>65\%$

Appliance service life: 1,000 cycles/treatments

Noise level at 50 cm*: 69.1 dB(A)

Service life of accessories: 1 year

* Data on noise levels were measured on a new device. Values may vary over time with use.

The information on the unit's capabilities provided by the manufacturer in accordance with EN 13544-1 may not apply to medicines provided in suspension form or those that are highly viscous.
 EN 60601-1 Medical electrical equipment – Part 1: general requirements for safety
 EN 60601-1-2 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance. Collateral standard: Electromagnetic compatibility - Requirements and tests.
 EN 13544-1 Respiratory therapy equipment – Part one: Nebulising systems and their components.

This section contains information specific to product compliance with the EN 60601-1-2 standard. The Chicco aerosol Super Soft is a medical electrical device that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned according to the electromagnetic information provided.

Mobile and portable RF communications equipment (mobile phones, transceivers, etc.) may affect the device.

The Chicco aerosol Super Soft is intended for use in the electromagnetic environment specified below. Customers or users of the product must ensure that it is used in the recommended environment.

Chicco nebulizers Super Soft [REF] 00009067000000 and [REF] 00009813000000 is intended for use in the electromagnetic environment specified below. The customer or the user of Chicco nebulizers [REF] 00009067000000 and [REF] 00009813000000 must make sure that it is used in such an environment.

Guidance and manufacturer's declaration - Electromagnetic emissions

Phenomenon	Professional healthcare facility environment ^{a)}	HOME HEALTHCARE ENVIRONMENT ^{a)}
Conducted and radiated RF EMISSIONS	^{a)}	CISPR 11 Group 1 Class B
Harmonic distortion	IEC 61000-3-2 ^{b)} Class A	
Voltage fluctuations and flickering	IEC 61000-3-3 ^{b)} COMPLIANT	

^{a)} The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The equipment is suitable for use in the mentioned environments when directly connected to the Public Mains Network.

^{b)} The test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard.

Guidance and manufacturer's declaration - Electromagnetic immunity - Enclosure port

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	Home healthcare environment
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8kV contact ± 2 kV, ±4kV, ±8 kV, ±15 kV air	

Radiated RF EM fields	IEC 61000-4-3	a)	10 V/m ^{b)} 80MHz - 2.7 GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	COMPLIANT NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the Chicco nebulizers Super Soft [REF] 00009067000000 (Frog) and [REF] 00009813000000 (Dog) can be requested from Artsana S.p.A. using the contact information provided in this manual. However, it is advisable to keep the electromechanical aerosol equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimise possible interference.	
RATED power frequency magnetic fields.	IEC 61000-4-8	30 A/m ^{c)} 50 Hz or 60 Hz	

a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied.

b) Before modulation is applied.

c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields.

Guidance and manufacturer's declaration - Electromagnetic immunity - Input a.c. port

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	Home healthcare environment
Electrical fast transients / Bursts.	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
Surges Line-to-line	IEC 61000-4-5	± 0.5 kV, ±1kV	
Surges Line-to-ground	IEC 61000-4-5	Class II equipment or supplied with Class II AC/DC adapter. Test not applicable	
Conducted disturbances induced by RF fields.	IEC 61000-4-6	3V ^{a)} 0.15 MHz - 80 MHz 6 V ^{a)} in ISM bands between 0.15 MHz and 80 MHz ^{b)} 80% AM at 1 kHz	
Voltage dips.	IEC 61000-4-11	0 % UT; 0.5 cycle At 0°,45°,90°,135°,180°, 225°, 270° and 315°.	
		0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0°	
Voltage interruptions ^{c)}	IEC 61000-4-11	0% UT; 250/300 cycles.	

^{a)} r.m.s., before modulation is applied.

^{b)} The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

^{c)} ME EQUIPMENT: the IMMUNITY TEST LEVEL is applied to the AC power input of the me equipment.

More information on distances to be kept between portable and mobile radio communications equipment and the Chicco Super Soft aerosol can be requested from Artsana S.p.A., through the contact information provided in this manual. It is nevertheless advisable to keep a safe distance (at least 1 m) between the aerosol therapy appliance and mobile phones or RF transmitters in order to reduce possible interference.

LEGEND OF SYMBOLS:

Warning! Follow the Instructions

I-O On - Off

Class II device

Hz Frequency

IP21

Device protected against penetration of solids and liquids.
Device protected against vertically dripping water.

SN Serial Number

REF Product identification code

Device with Type BF applied part

Manufacturer

CE 0068

Complies with Directive MDD
93/42/EEC + 2007/47/EC

~ Alternating current

Device date of production: the second pair of digits in the lot number identifies the year, the third pair identifies the month. (**LOT** 1016100001; 16 = 2016; 10 = October)

CONDITIONS FOR STORAGE AND TRANSPORTATION:

Relative humidity: 15% / 85% RH

Temperature: -25°C / 70°C

Atmospheric pressure from 700hPa to 1060hPa

ENVIRONMENTAL CONDITIONS DURING USE:

Relative humidity: 15% / 85% RH

Temperature: 5°C / 40°C

Atmospheric pressure from 700hPa to 1060hPa

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions. The warranty shall not therefore apply in the case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national laws in the country of purchase, where appropriate.

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Artsana S.p.A - Via Saldarini Catelli, 1