EASY BREATH GO+

AEROSOL MESH CON MEMBRANA SOSTITUIBILE
MESH NEBULISER WITH REPLACEABLE MEMBRANE
NÉBULISEUR MESH AVEC MEMBRANE REMPLAÇABLE
MESH-NEBULISATOR MIT AUSTAUSCHBARER MEMBRAN
NEBULIZADOR MESH CON MEMBRANA REEMPLAZABLE
NEBULIZADOR MESH COM MEMBRANA SUBSTITUÍVEL
NEBULIZATOR MESH Z WYMIENNAJ MEMBRANĄ
MESH-NEBULISATOR MET VERVANGBARE MEMBRAAN

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USER MANUAL

Dear customer, thank you for choosing Chicco Easy Breath on the Go+, a mesh nebuliser designed and manufactured according to the most up-to-date technologies. The nebuliser is a class IIa medical device according to European Regulation for medical devices EU Reg 2017/745. Before using the device, we highly recommend to read this manual, to make sure you have correctly understood how to use the nebulizer.

IMPORTANT: Read the instructions carefully before use.

Mesh unit (2b): please be sure that it is well assembled before starting using the device.

Do not boil the mesh membrane.

Do not use brush, cotton buds, needles or other things to clean or dry the mesh unit, this could cause irreparable damage.

Clean and disinfect following the instructions in this manual.

Device not suitable for nebulizing essential oils, phyto-therapeutics, hyaluronic acid, hypertonic solution, thermal waters naturally rich in salts.

Remove the batteries if the product remains unused for a long time.

COMPONENTS

1. Cover of medication container with silicone o-ring
2. Medication container, composed of medication tank (2a), mesh unit (2b) and connector (2c)
3. Body of nebuliser
4. On/off button
5. Power LED indicator
6. Maintenance LED indicator
7. Button for unlocking the medication container
8. USB-C port
9. Battery cover
10. Child mask
11. Adult mask
12. Mouthpiece
13. USB cable
14. Batteries
15. Spare mesh unit

INTENDED PURPOSE:

The device is an aerosoltherapy system suitable for domestic use.

The device is designed to drive the piezoelectric ceramics mesh nebulizer module to produce medical aerosol for lung respiratory disorders.

PATIENT POPULATION:

The nebuliser can be used on person over 2 years of age under supervision; it can be used for self- treatment by persons over 12 years of age.

INTENDED USERS:

The user is an adult person without specific knowledge or professional ability. The patient is the intended operator except in the case of children and patients that requires special assistance.

INDICATIONS:

Lung acute or chronic diseases of respiratory disorders organs, or inflammation of the upper respiratory system.

CONTRAINDICATIONS:

Device is not indicated to use with quick-relief medications during life-threatening asthma attacks.

No contraindication exists to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

WARNINGS:

Please report any serious incident that has occurred in relation to the device injury or adverse event to the competent local authority and to the Manufacturer or to the European Authorized Representative (EU REP) Vigilance contact points: https://health.ec.europa.eu/medical-devicessector/new-regulations/contacts_en.

- Use the device only as described in this manual and therefore as an aerosol therapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- The device is not intended for health care in professional facilities.
- If children or people with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device.

- Do not use the device on people who are ventilated and/or unconscious.
- Before use, ensure that there is no visible damage to the device or components
- Do not use the device, if you spot damage or you notice something unusual, if the device does not work properly or if it's making an abnormal sound. If you feel unwell or experience pain, stop using it immediately.
- In case of failure and malfunction, read the "Troubleshooting - Malfunctions and Actions to take" in the user manual. Do not handle or open the device
- Comply with the safety regulations concerning electrical devices and in particular:
 - use only original applied parts;
 - never submerge the unit in water;
 - do not use the device while taking bath or shower
 - never wet the device, it is not protected against water penetration;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
- Before plugging into the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated or has fully charged.
- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to people, animals or things, for which the manufacturer cannot be held responsible.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply or remove batteries.
- Some parts of the unit are so small that they may be swallowed by children; keep the equipment out from children's reach
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
 - use this device only with medicines prescribed by your doctor;
 - only administer the treatment using the applied part recommended by your doctor depending on the pathology;
- Check the medicine package leaflet for possible contraindications for use with common aerosol therapy systems
- The device is not a substitute for medical consultation and treatment. Always consult your doctor first in case of pain or illness.
- Observe the general hygiene measures when using the nebuliser.
- To avoid strangulation and entanglement, keep cable out of reach of young children.
- Do not position the equipment so that it is difficult to operate the disconnection device.
- Nebuliser and applied parts are single-patient reuse. Device is multi-patient reuse.
- Observe the general hygiene measures when using the nebuliser.
- Changes to the device and the components are not permitted.
- Use only applied parts recommended by the manufacturer, use of different components may impair the effectiveness of treatment and damage the device.
- Under no circumstances should you open or repair the device yourself. There is a risk of electric shock when opening the device. The device is not disconnected from the power supply unless the adapter is unplugged (and the USB-C cable has no other power connection).
- Never use the device near flammable gases, explosive gas, oxygen or nitrogen oxide.
- Do not operate the device near electromagnetic transmitters.
- Protect the device from impacts.
- Should you require an adapter or extension lead, it must meet the applicable safety requirements. The power limit and the maximum output specified on the adapter must not be exceeded.

HOW TO USE THE DEVICE

The device must be checked before each use in order to identify any malfunctions and/or damage caused during transportation and/or storage. During inhalation, sit upright and relaxed to avoid compressing the airways and impairing treatment effectiveness. The applied parts must be used only with a single patient, it is not recommended to use them with several patients.

After unpacking the device, check it for visible damage or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check if the applied parts are intact.

Before using the device, proceed with the cleaning operations as described in the chapter "CLEANING AND MAINTENANCE".

1. Preparing the nebulizer:

Assemble Mesh unit to the medication container: pair mesh unit (2b) with medication tank (2a) by inserting the mesh unit with electrical contacts downward (C1), then close by applying the connector (2c) and rotate clockwise (C1), until the arrow indications of connector and tank are aligned (C3).

Check that mesh unit is correctly inserted in medication cup, by paying attention that mesh metal component is visible (C2A) from outside, otherwise liquids may leak from the device.

BATTERY MODE: open the battery cover, insert the batteries into device by paying attention to polarity and close again as shown in pictures B1 and B2.

USB – MODE: connect USB-C cable (13) to the device port (8) and to the USB-A power supply (not included in delivery. The external power supply must comply with the European standard EN 60601-1 with protection class 2 and must have the rating specified in section "Technical specifications").

Filling the medication container

- Open the cover of the medication container (D) and fill the medication container (E) as needed. The maximum recommended filling quantity is 10 mL. Avoid overfilling, otherwise liquids may leak from the device.
- Use medication only as instructed by your doctor and ask about the appropriate inhalation period and quantity for you.
- If the prescribed quantity of medication is less than 0.5 mL, it must be topped up with isotonic saline solution. Dilution is also necessary with viscous medications (viscosity < 3). Observe the instructions of your doctor.
- Close the cover of the medication container (F) until the markings on cover and medication container are aligned.

2. Fix applied parts and start therapy

- Securely connect mask or mouthpiece to the connector (G). If you are using the adult mask or child mask, you may attach one of the straps included. Ensure that you securely fix the strap on the mask with a knot.
- Now move the device to your mouth and put your lips firmly around the mouthpiece. When using a mask, place it over nose and mouth, without covering the eyes.
- Start the device by pressing the On/Off button (4).
- If the spray mist is flowing out of the device and the power LED (5) is lit up the device is operating. Refer to "LED display color" section for details about LED color in different operation modes.
- Keep the device as vertical as possible. However, a small angle does not affect usage as the container is leak-proof. Make sure that the medication is in contact with the mesh for proper functioning of the device during use. When liquid is about to be run out user is suggested to slightly incline the device toward himself or herself.

3. End of nebulization

- If you have finished or want to stop the treatment, you can switch off the device manually using the On/Off button (4)
- The device has an automatic switch-off function, which is activated:
 - After a nebulisation period of 20 minutes.
 - When there is only a very small amount of medication or liquid left or when the substance to be nebulized is no longer in contact with the mesh. The power LED (5) flashes quickly and the device then switches off.
- Do not use the device if the medication container is empty. The device detects when there is no medication in the medication container and will switch-off automatically.
- For technical reasons, a small amount of medication will remain in the medication container. Do not reuse it.

LED display color

Power LED (5) – Maintenance LED (6)

LED DISPLAY COLOR	MEANING
Power LED Green	Device in USB MODE
Power LED Blue	Static ON: device in BATTERY MODE Fast Blinking for 3 seconds: absence of liquid inside the liquid container
Power LED Orange	Fast Blinking: device is running in Auto-cleaning mode Static ON: low charge battery
Power LED Green-Blue-Orange	Blinking for 3 seconds: medication container is missing
Maintenance LED Orange	Static ON: reminder signal that the mesh unit must be changed Fast blinking: informatory signal that Auto-cleaning is needed

CLEANING AND MAINTAINING THE DEVICE

Before the first use and after each treatment

 Follow carefully the cleaning and disinfection instructions as they are very important for the device performances, the therapy success and to avoid health risks.

- Do not clean the mesh mechanically using a brush or similar device, this could cause irreparable damage.
- Please consult your doctor about the additional requirements for hygiene preparation (hand care, handling of medication/inhalation solutions) for high-risk groups (e.g. patients with cystic fibrosis).

Limitations on processing

The maximum time of cleaning and disinfection of medication container (2a), connector (2c), mouthpiece (12) and mask(s) (10-11) is 366 times.

To maintain proper device functioning:

- We recommend that you replace the medication container, mouthpiece and mask(s) after maximum number of cleaning cycles, depending on the usage.

Mesh unit replacement

We recommend to replace mesh unit (2b) after 1000 operating minutes. Device provides reminder signal for

mesh unit replacement. To replace the mesh unit, follow disassembling and reassembling instructions provided in "Preparation before cleaning" section. Before reassembling the medication container, timer can be reset by pressing on/off button for 5 seconds from device off condition.

INSTRUCTIONS	
Preparation before cleaning	<p>Disassembling the device</p> <ul style="list-style-type: none"> - Press the release button on the rear of the device and vertically pull off the medication container (H). - Remove the mouthpiece or the mask from the medication container - Open the cover of the medication container (D) - Separate connector (2c), mesh unit (2b) and medication tank (2a) (C1) by rotating at first the connector counterclockwise and then slightly tilting the tank forward to pull out the mesh unit. <p>Note: To avoid microbial growth and residual drug drying, clean and disinfect the device immediately after each use.</p>
Cleaning	<p>After each use, the disassembled medication container, connector, mesh unit, mouthpiece or mask, must be cleaned by immersion in a clean container with soft warm water at a max. temperature of 40°C for approx. 5 minutes. Scrub all the surfaces of all the components (except the mesh) with a clean and small brush for at least 8 times. Afterwards rinse all parts (except the mesh unit) thoroughly in running tap water for at least 30 seconds</p> <ul style="list-style-type: none"> - Do not clean the mesh unit mechanically using a brush or similar device - Do not wash the mesh unit under running water. - Do not remove silicone O-ring from cover (1) <p>Reassembly is carried out in reverse order, as indicated in paragraph How to use the device.</p>
Disinfection	<p>Place the dissembled medication container, the mesh unit, the connector, mouthpiece and masks in a clean container. Cover everything with ethyl alcohol (70-75%) and leave it for 10 minutes. Finally, clean all parts again with sterile water (soft water) and allow to air dry in a clean place.</p>
Disinfection (alternative not for mesh unit)	<p>Alternatively, medication tank (2a), connector (2c), cover (1), mouthpiece (12) and masks (10-11) can be placed in boiling water for 5 minutes for disinfection.</p> <p>Note: The mesh must not be placed in boiling water.</p>
Drying	<ul style="list-style-type: none"> - Gently shake the parts a few times so that any water is removed from the tiny holes in the mesh. - Place the individual parts on a clean, dry and absorbent surface and leave them to dry completely (for at least 4 hours). <p>Note: do not use microwave for drying.</p>
Inspection	<p>Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts.</p>
Packaging	<p>Pack dry parts in a clean and sealed container when not in use. Do not pack wet or damp parts.</p>
Storage	<p>Storage conditions refer to "TECHNICAL SPECIFICATIONS"</p> <p>Note: re-clean and re-disinfect the parts if they are stored more than one day.</p>
Transportation	<p>After cleaning and disinfection, always transport the parts within clean and sealed container. Contamination can be avoided by good hand washing and not touching the inside sections of the parts when take out and re-assemble the parts for use.</p>

Auto-cleaning function for medication container

Auto-cleaning is recommended periodically. Device provides reminder for auto-cleaning after 300 operating minutes. Timer is reset when auto-cleaning function is activated.

- Open the cover of the medication container (D) and pour out the remaining fluid.
- Pour a sufficient amount (3-6 ml) of tap water (soft water) into the medication container (E). Then close the cover (F). Shake the container slightly so that the remaining fluid mixes properly with the cleaning fluid. Let the device run for one minute and then switch it off.
- Open the cover of the medication container (D) and pour out all the liquid.
- Pour 1.5 to 2.0 ml of tap water or distilled water mixed with white vinegar (1/4 vinegar and 3/4 distilled water) into the connector (2c) of the medication container placed in horizontal position (I). To avoid leakage through the lateral openings, please do not overfill and hold the body firmly.
- From OFF status press and hold the On/off button for 3 seconds: the power LED (5) flashes blinking orange and the product switches to auto-cleaning mode. The cleaning liquid is now drawn out from the connector (2c) into the medication tank (2a) through mesh unit (2b) in reverse mode. During this process, it is possible for droplets of this liquid to

leak out of the medication container. Therefore, place some protection under the nebulizer if necessary. Auto-cleaning will automatically run for 3 minutes.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

TROUBLESHOOTING

PROBLEM/QUESTION	POSSIBLE CAUSE/REMEDY
The nebuliser produces no or too little aerosol.	<ul style="list-style-type: none"> • Too little medication in the medication container. • The nebuliser is not being held vertically. • Medication liquid has been poured in that is not suitable for nebulisation (e.g. too viscous. The viscosity must not be greater than 3. Ask your pharmacist about the viscosity.) The medication liquid should have been prescribed by a doctor. • The battery is flat, change the the batteries or use with USB cable. • The mesh unit might be damaged/worn. It is recommended to replace it. • Particles on the mesh are hindering output. In this case, you should pause the inhalation process and remove the mask or mouthpiece. Then clean and disinfect the medication container (see "Cleaning and maintaining" section). • The mesh is worn. Replace the mesh unit.
What medications are suitable for inhaling?	<p>Only your doctor can decide which medication should be used for treating your condition. Please consult your doctor in relation to this matter.</p> <p>With the device, you can nebulise medications with a viscosity of less than 3. Medications containing oil (particularly essential oils) must not be used.</p>
There is inhalation solution residue in the medication container	This is normal and due to technical reasons. Stop the inhalation as soon as the nebuliser makes a significantly different noise or the device switches off automatically due to a lack of liquid.

TECHNICAL SPECIFICATIONS

ITEM	SPECIFICATION
Type	GCEC03
Model No.	Easy Breath Go+
Supply	2x1.5V AA batteries/Input: 5V DC(USB),500mA
Filling volume	min. 0.5 ml, max. 10 ml
Sound pressure	approx. 42 dB(A)
Vibration frequency	108 kHz \pm 10%
Aerosol values according to EN ISO 27427 based on adult breathing patterns with sodium fluoride (NaF):	<ul style="list-style-type: none"> • Aerosol output: 0.28 ml • Aerosol output rate: 0.05 ml/min • Residual volume: 0.63 ml • Particle size (MMAD): 3.89 \pm 0.35 μm • GSD (geometrical standard deviation): 11.81 \pm 0.05 • RF (respirable fraction < 5 μm): 63.0 \pm 5.86% • Large particle range (> 5 μm): 37.0 \pm 5.86% • Medium particle range (2 to 5 μm): 48.3 \pm 1.87% • Small particle range (< 2 μm): 14.7 \pm 4.07%
Nebulisation rate (average value) in USB mode:	0,4 ml/min
Nebulisation rate (acceptable range) in USB mode	0,25–0,7 ml/min
Expected service life of the device	device: 500 hours, mesh unit: 1000 min
Equipment conditions of use	20 min ON / 40 min OFF
Operation condition	10°–40°C, 15–90% RH, 700hPa - 1060hPa
Storage condition	-20°–60°C, 15–90% RH, 700hPa - 1060hPa

Transport condition	-20°~60°C, 15-90% RH, 700hPa - 1060hPa
---------------------	--

Device complies with European Regulation on Medical Devices EU MDR 2017/745.
Class II device referred to protection against electric shock.

Important for the safety of the device and its user; use only power supplies as specified (5V > 500mA output with USB-C connector) and ensure their certification in accordance with EN 60601-1 (with protection class 2) and EN 60601-1-2. Note: Technical specifications subject to change without notice.

NOTES ON ELECTROMAGNETIC COMPATIBILITY

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be restricted if electromagnetic disturbance occurs. This could result in issues such as the failure of the device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of an atomiser set other than that specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment, such as antenna cables or external antennas) at least 30 cm away from all device parts, including all cables included in delivery.
- Failure to comply with the above could impair the performance of the device.

MAINTENANCE AND REPAIRS

In case of failure, contact Customer Service. **Never open or try to repair the device.**

The device has no user-serviceable parts and does not need internal maintenance or lubrication.

DISPOSAL

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

THIS PRODUCT COMPLIES WITH EU REGULATION 2023/1542

The crossed bin symbol on the battery or on the packaging means that the battery at the end of its life should be disposed separately from the municipal waste stream via designated collection facilities appointed by government or local authorities or otherwise returned to the dealer when purchasing new equivalent rechargeable or non-rechargeable battery. The chemical symbols Hg, Cd, Pb printed under the cross-out wheeled bin symbol indicates the type of substance contained in the battery. Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for taking the battery to a suitable waste disposal site at the end of product's life for processing and recycling. If dead batteries are collected properly as separate waste they can be recycled. Processed and ecologically disposed; this avoids a negative impact on both environment and human health and contributes towards the recycling of the substance contained in the batteries. Improper disposal of the product by the user may harm the environment and human health. For further information regarding waste disposal contact your local waste disposal agency or the shop where you bought the appliance.

For full product disposal description please refer to

<https://www.chicco.com/medical-device.html> and <https://www.chicco.com/recycle.html>

WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions. The warranty shall not therefore apply in case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national laws in the country of purchase, where appropriate.

DEVICE SYMBOLS - MANUFACTURER'S DECLARATION AND GUIDELINES - ELECTROMAGNETIC EMISSIONS: SEE APPENDIX.

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DEVICE SYMBOLS

	Reference Number		Distributor
	ON /OFF		Authorized representative
	Direct current		Manufacturer
	Class II		Date of manufacture (year, month and day)
	Type BF		Serial Number
	Device Complys with European Regulation on Medical Devices EU MDR 2017/745		Caution
	Medical device		Humidity limitation
	Follow instruction for use		Temperature limitation
	Importer		Ambient pressure limitation
	Unique device identifier		Type Number
	Batch code		Single patient multiple use (for applied parts only)
	Protection against solid foreign objects and against harmful effects due to ingress of water		Authorized representative

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 2\text{kV} \pm 4\text{kV} \pm 6\text{kV} \pm 8\text{kV}$ $\pm 15\text{kV}$	$\pm 8\text{ kV}$ contact $\pm 2\text{kV} \pm 4\text{kV} \pm 8\text{kV} \pm 15\text{kV}$ air	The floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	$\pm 1\text{kV}, \pm 2\text{kV}, 100\text{ kHz}$ repetition frequency	Power supply lines: $\pm 2\text{ kV}$	The main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0.5\text{kV}, \pm 1\text{kV}$ (Line to line) $\pm 0.5\text{kV}, \pm 1\text{kV}, \pm 2\text{kV}$ (Line to Ground)	Line to line: $\pm 0.5\text{kV}, \pm 1\text{kV}$ Line to Ground: $\pm 0.5\text{kV},$ $\pm 1\text{kV}, \pm 2\text{kV}$	The main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U^T for 0.5 cycle 0% U^T for 1 cycle 70% U^T for 25/30 cycles 0% U^T for 250/300 cycles	0% U^T for 0.5 cycle 0% U^T for 1 cycle 70% U^T for 25/30 cycles 0% U^T for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) Magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U^T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Smart including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. Recommended separation distances $d = 1,2 \sqrt{P}$ $d = 2 \sqrt{P}$

Radiated RF IEC61000-4-3	10V/m	10V/m	<p>80 MHz a 800MHz $d=1,2 \sqrt{P}$ Da 800 MHz a 2.7 GHz $d=2,3 \sqrt{P}$</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a) Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the device

The Device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d=1.2\sqrt{P}$	150 kHz to 80 MHz (in ISM and amateur radio bands) $d=2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.7GHz $d=2.3\sqrt{P}$
0,01	0.12	0.2	0.12	0.23
0,1	0.38	0.632	0.38	0.73
1	1.2	2	1.2	2.3
10	3.8	6.32	3.8	7.3
100	12	20	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

	Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{b)}	Modulation ^{c)}	Modulation ^{d)} (W)	Distance (m)	Immunity Test Level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation ^{e)} 18 Hz	1,8	0,3	27
	450	380 – 390	GPRS 460, FR5 460	FM ^{f)} ± 5 kHz	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation ^{e)} 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{e)} 18 Hz	2	0,3	28
	870						
	930						
	1 720	1 700 – 1 990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation ^{e)} 217 Hz	2	0,3	28
	1 845						
1 970							
2 450	2 400 – 2 570	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation ^{e)} 217 Hz	2	0,3	28
	5 240						
	5 500	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation ^{e)} 217 Hz	0,2	0,3	9
	5 785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50% duty cycle square wave signal.

^{c)} As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The Manufacturer should consider reducing the minimum separation distance, based on Risk Management, and using higher Immunity Test Levels that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher Immunity Test Levels shall be calculated using the following equation:

$$E = 6/d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the Immunity Test Level in V/m

NOTE

A series of horizontal dotted lines for writing notes, spanning the width of the page.

NOTE

A series of 18 horizontal dotted lines for writing notes.

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REF 00011481000100

Type: GCEC03

Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road,
European Industrial Zone, Xiaolan Town, 528415 Zhongshan City,
Guangdong Province, PEOPLE'S REPUBLIC OF CHINA,
Phone: +86 760 22589901, email: service@globalcare.com.cn

Donawa Lifescience,
Piazza Albania, 10, 00153 Rome, Italy
<https://www.donawa.com/wli/main/contatti.index>

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Artsana S.p.A