# EASY BREATH BABY

Aerosol a pistone

Piston Nebulizer

REF 00011482000000


FIG. A



12

13

14

## **II** ISTRUZIONI D'USO

Gentile cliente,

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# EN USER MANUAL

Dear customer, thank you for choosing this Chicco nebuliser piston, an nebuliser therapy unit designed and manufactured according to the most up-to-date technologies. The nebuliser is a class IIa medical device according to European Regulation for medical devices EU Reg 2017/745. The device can be used together with Easy Breath Nasal Shower [REF](#) 00012247000000, ideal for hydrate and wash nasal cavities. Additional Easy Breath Kit [REF](#) 00012248000000 may be purchased separately. Before using the unit, we highly recommend that you read the shortlist of warnings contained in this manual, to make sure you have correctly understood how to use it.

## IMPORTANT:

**Read the instructions carefully before use.**

## COMPONENTS

1. Piston compressor
2. Power lead
3. ON/OFF Switch
4. Air filter compartment
5. Ampoule (nebulization cup)
6. Air tube
7. Mouthpiece
8. Nasal fork
9. Adult mask
10. Paediatric mask
11. Air filters
12. Removable accessory compartment
13. Mask-ampoule connection
14. Reassuringly nebuliser cap

## INTENDED PURPOSE:

The device is an nebuliser therapy system suitable for domestic use.

The device is designed for the production of compressed air to operate a small volume nebulizer for the production of medical nebuliser for lung respiratory disorders.

## PATIENT POPULATION:

The device is intended for use with children from 2 years old, adolescent and adult patients.

## INTENDED USERS:

The user is an adult person without specific knowledge or professional ability. The patient is the intended operator except in case of children and patients that requires special assistance. If children or persons with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device.

## INDICATIONS:

Lung acute or chronic diseases of respiratory disorders organs, or inflammation of the respiratory system.

## CONTRAINDICATIONS:

The device is not indicated to use with quick-relief medications during asthma attack. No contraindications exist to the administration of nebulisers by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

**WARNINGS:**

Please report any serious incident that has occurred in relation to the device injury or adverse event to the local competent authority and to the Manufacturer or to the European Authorized Representative (EC REP) Vigilance contact points: [https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts\\_en](https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en)

- Use the device only as described in this manual and therefore as an nebuliser therapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- Keep the user manual for future reference.
- Do not operate the unit in presence of any anaesthetic mixture inflammable with oxygen or nitrogen protoxide.
- The correct functioning of the equipment can be affected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it to a different power socket.
- In case of failure and/or malfunction, read the "Malfunctions and Actions to take" in the user manual. Do not handle or open the compressor housing.
- For repair operations address only to a technical service center authorized by the manufacturer and require the use of original spare parts. The non-observation of the above mentioned indications can compromise the device safety.
- Comply with the safety regulations concerning the electrical devices and in particular:
  - use only original accessories and components;
  - never submerge the unit in water;
  - never wet the device, it is not protected against water penetration;
  - never touch the unit with wet or moist hands;
  - do not leave the unit exposed to the weather elements;
  - place the unit on a stable and horizontal surface during its operation;
  - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
  - do not pull the power cord or the device itself to unplug it from the power socket;
  - the power plug is the separation element from the grid power; keep the plug accessible when the device is in use.
- Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- Use of adapters, simple or multiple, and/or extension cables is not recommended. If their use is indispensable, it is necessary to use types complying with safety regulations, paying however attention that they do not exceed the maximum power limits, indicated on adapters and extension cables.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated.
- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or things, for which the manufacturer cannot be held responsible.
- The power cord of this device cannot be replaced by the user. In case of a power cord damage, address to a technical service center authorized by the manufacturer for its replacement.

- The power supply cord should always be fully unwound in order to prevent dangerous overheating.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply.
- Keep the equipment out from children's reach: some parts of the unit are so small that they may be swallowed by children and the power cord poses risk of strangling.
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
  - use this device only with medicines prescribed by your doctor;
  - make the treatment using only the accessory recommended by your doctor depending on the pathology;
  - use the nosepiece accessory only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose, but only bringing them as close as possible.
- Check in the medicine package leaflet for possible contraindications for use with common nebuliser therapy systems.
- To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
- Do not position the equipment so that it is difficult to operate the disconnection device.
- Nebuliser and accessories are single patient reuse. Device is multipatient reuse
- Nebulizing system is not suitable for use in anaesthetic breathing system or a ventilator breathing system. The patient cannot undergo MRI scan while using this device.
- Device is designed to nebulize solution and suspension liquids
- Do not use the device, if you spot damage or you notice something unusual.
- This device consists of delicate components and must be handled with care. Observe the storage and operating conditions in the chapter "Technical Data"
- The use of this device is not a replacement for visiting the doctor. In case of no health improvement after treatment consult again your doctor.

## HOW TO USE THE DEVICE

The device must be checked before each use in order to identify any malfunctions and/or damage caused during transportation and/or storage. During inhalation, sit upright and relaxed to avoid compressing the airways and impairing treatment effectiveness. The accessories must be used only with a single patient, it is not recommended to use them with several patients. After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check if the accessories are intact. Before using the device, proceed with the cleaning operations as described in the chapter «CLEANING AND MAINTENANCE». For an even quieter therapy, remove the detachable accessory holder container and the cap during nebulization.

1. Open the nebuliser cup by turning counter clockwise the top (Pic. B)
2. Make sure that the medicine conduction cone is properly fitted on the air conduction cone inside the nebuliser cup (Pic. C).
3. Pour the medicine and/or saline solution in the nebuliser cup, according to the doses recommended by your physician, without exceeding the maximum level of the nebuliser cup
4. Close the nebuliser cup by turning clockwise the two parts, paying attention that they are thoroughly sealed (Pic. E).
5. Connect one end of the air tube to the nebuliser (Pic. F) and the other end to the air outlet

- (Pic. G) on the device.
6. Connect the accessory needed for the nebuliser therapy (mask, mouthpiece or nasal fork) (Pic.H)
  7. Plug the device into the wall socket, making me that the main supply corresponds to the electrical rating of the device.
  8. To start the treatment, set the I/O switch into the "I" position. Inhale the nebuliser solution using the prescribed accessory.
  9. When the treatment has been completed, switch the unit off by setting the I/O into the "O" position and disconnect the plug from the wall socket.
  10. Remove the remaining medication from the nebuliser cup and clean the device as described in the section «Cleaning and Disinfecting».
  11. The device requires no calibration.
  12. Do not block the ventilation hole located at the lower part of the device. It may case the device heat up or break or case injury.
  13. Using the product continuously may heat it up.

This device is for intermittent use: 20 min ON / 40 min OFF.

Do not use the product continuously for 20 minutes or longer.

Using the device before 40 minutes have elapsed may result in motor overheating and cause the safety cut-off devices to activate.

To prevent overheating and damage to the compressor, switch off the appliance for at least 40 minutes after every 20 minutes of uninterrupted use.

## CLEANING AND DISINFECTION

Device(s) : nebuliser for aerosoltherapy and accessories	
WARNINGS	<p><b>Clean and disinfect before the first use and after each treatment.</b></p> <p>Follow carefully the cleaning and disinfection instructions of the accessories as they are very important for the device performances and the therapy success. <b>Use original accessories only. Do not clean or disinfect the air tube. The accessories cannot be cleaned and disinfected by automated method.</b></p>
Limitations on processing	<p>The nebuliser must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebuliser nozzle is obstructed by dry medicine, dust, etc. We recommend to replace the nebuliser after a period between 6 months and 1 year, depending on the usage. The maximum times of cleaning and disinfection of the nebuliser are 360 times. The maximum times of cleaning and disinfection of the mask, the nosepiece and the mouthpiece are 360 times.</p>

<b>INSTRUCTIONS</b>	
Preparation before cleaning	Detach the air tube from the nebuliser. Detach the mouthpiece or nosepiece or mask from the nebuliser. Make sure that all residual volume is removed from the nebuliser. Disassemble the nebuliser cup by turning counterclockwise the top and remove the medicine conduction cone. Note: To avoid microbial growth and residual drug drying, clean and disinfect the accessories immediately after each use
Cleaning	Briefly rinse all parts in running tap water for at least 10 seconds beforehand. Mix little dish soap with warm water in a clean container. water i.e. Immerse the components of disassembled nebuliser cup, the mask, the mouthpiece and the nasal fork into the water with soap for about 5 min. And then scrub all the surface of all the components with clean small brush for at least 8 times. Afterwards rinse all parts thoroughly in running tap water at least 30 seconds, to fully remove any possible dish washing liquid residue.
Disinfection	After cleaning, disinfect all of the disassembled parts (only parts that have been cleaned can be disinfected effectively). Boil the disassembled nebuliser cup, the mouthpiece and the nasal fork for 5 minutes in boiling water. DO NOT BOIL THE MASK. Immerse the mask in sodium hypochlorite (NaOCl) solution for 15 minutes, (i.e. a solution made from the disinfecting agent). Afterwards immerse the mask with sterile water for 3 min, and then rinse all the surface of the mask with other clean sterile water twice, to fully remove any possible residue of the disinfectant solution.
Drying	Reassemble the nebuliser components and connect it to the air-outlet, switch the device on and let it work for 10-15 minutes. Let all parts dry completely before reassembling and reuse to avoid risk of germ grow. Care needs to be taken not to contaminate the parts after they have been cleaned and disinfected. Contamination can be avoided by good hand washing and not touching the inside sections of the device when laying them out to dry or when reassembling.
Inspection	Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts.
Packaging	Pack dry parts in a clean and sealed container when not in use. Do NOT pack wet or damp parts.
Storage	Storage conditions refer to "TECHNICAL SPECIFICATIONS". Note: re-clean and re-disinfect the parts if they are stored more than one day.

Transportation	After cleaning and disinfection, always transport the parts with clean and sealed container. Contamination can be avoided by good hand washing and not touching the inside sections of the parts when take out and re-assemble the parts for use.
----------------	---

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the user to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### **FILTER REPLACEMENT**

Periodically check the condition of the filter. The filter was inserted to protect the compressor. Correct filter maintenance prolongs the life of the nebuliser. The frequency with which filter should be replaced depends on the conditions in which the device is used. In normal usage conditions, the air filter must be replaced approximately after 200 working hours or once a year. We recommend to periodically check the air filter (10-12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Remove the filter and replace it with a new one. Do not attempt to clean the filter to reuse it. The air filter shall not be serviced or maintained while in use with a patient.

The filter is located on the bottom of the device.

#### **TO REPLACE THE FILTER:**

- Remove the filter cover cap turning it anti-clockwise;
- Remove the filter to be replaced from its housing (using a small tool if necessary);
- Carefully fit the new filter in its housing;
- Replace the filter cover cap and fasten tightly by turning it in a clockwise direction;
- **WARNING!** Never leave the removed filter cover cap and/or filters unattended: These are small parts posing the risk of suffocation if swallowed by children.

**Only use original filters. Do not use the device without the filter**

### **MALFUNCTIONS AND ACTIONS TO TAKE**

The device cannot be switched on:

- Ensure the power lead is correctly plugged into the socket.
- Ensure the ON/OFF switch is in the position «I».
- Make sure that the device has been operating within operating limits indicated in this manual (20 min On / 40 min Off).

The nebuliser functions poorly or not at all:

- Ensure the air tube (6) is correctly connected at both ends.
- Ensure the air tube is not squashed, bent, dirty or blocked. If necessary, replace with a new one.

For other nebulization problems:

TYPE OF DEFECT	CAUSE	REMEDY
Low nebulization	Clogged nebuliser cup	Remove and detach all the nebuliser cup components and clean them as illustrated in this instruction manual.
Low nebulization	Clogged nebuliser cup	If the issue persists after cleaning, replace the nebuliser cup.
No nebulisation	Nebuliser cup components not assembled correctly	Disassemble and reassemble the nebuliser cup; if the issue persists, replace the nebuliser cup.
No nebulisation	The compressor does not work	Check that the power cord plug is disconnected from the mains and that the on/ off switch is on. if the issue persists contact the customer service.
Slow nebulisation	Drug too thick	Dilute the drug as instructed by the manufacturer.

In the event that the equipment does not resume correctly, please contact **CUSTOMER SERVICE**.

## MAINTENANCE AND REPAIRS

In case of failure, contact **CUSTOMER SERVICE**. **Never open the device**. The device has no user-serviceable parts within and does not need internal maintenance or lubrication.

## TECHNICAL SPECIFICATIONS

Item	Specification	
Max Pressure	200-320	kPa
Operating Pressure (with calibrated nozzle, $\Phi$ 0,50)	80 – 155	kPa
Protection type IP	IP 21	--
Sound noise level (Distance: 1m, Front direction) max level	$\leq 58$	dB(A)
Equipment conditions of use	20 ON / 40 OFF	min
Dimensions	177.1x177.1x176.8	mm
Weight	1.7 $\pm$ 0.2kg	kg
Service life	1000	Hours
Power supply	230V AC 50Hz	--
Neb Rate (Nominal value)	0.35	ml/min
Neb Rate (acceptable range)	0.23~0.60	ml/min
Nebuliser kit NEBULISER PERFORMANCES ACCORDING TO EN ISO 27427:2019 based on adult ventilatory pattern with Sodium Fluoride (NaF)	Jet Baby	--

Output	0.11	ml
Output rate	0.07	ml/min
Percentage of fill volume emitted per min	3.5%	-
Residual Volume	1.68	ml
Particle size (MMAD)	3.07	µm
GSD (geometric standard deviation)	2.03	-
RF (respirable fraction <5µm)	75%	-
Large particle range (>5µm)	25%	-
Middle particle range (2 to 5µm)	47.8%	-
Small particle range (<2µm)	27.2%	-
Min & Max filling volume	2 - 8	ml
Operating Temperature	10 ~ 35	°C
Operating Humidity	10 ~ 95	% RH
Operating Pressure	700 ~ 1060	hPa
Storage Temperature	-25 ~ 70	°C
Storage Humidity	10 ~ 95	% RH
Storage Pressure	700 ~ 1060	hPa

- Class II device referred to protection against electric shock.
- Nebuliser, mouthpiece, masks, and fork are type BF applied parts Note: Technical specifications subject to change without notice.

### NEBULISER CHARACTERISTICS IN ACCORDANCE WITH REGULATION EN ISO27427 Standards applied:

Electric safety standards CEI EN 60601-1 Electromagnetic compatibility according to CEI EN 60601-1-2 The device is a Class IIa medical device.

### DISPOSAL

ITEM	COMPONENT	DISPOSAL DESCRIPTION	PHOTO
1	Device	The component mainly includes plastic and electronic components. All comply with RoHS and REACH, and all could be safety disposal. This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products.	

2	Accessories (tube, nebuliser bottle, nosepiece, mouthpiece, mask)	The component are plastic. All comply with ROHS and REACH. Before disposal, all the accessories should be cleaned according to manual instructions and then disinfected.	
---	---	--	---

## DEVICE SYMBOLS

	Catalogue Number		Distributor
	ON /OFF		Authorized representative
	Alternating current		Manufacturer
	Class II		Date of manufacturing (year, month and day)
	Type BF		Serial Number (LOT/serial number)
	Device Complies with European Regulation on Medical Devices EU MDR 2017/745		Caution
	Medical device		Humidity limitation
	Follow instruction for use		Temperature limitation
	Importer		Ambient pressure limitation
	Unique device identifier		Type Number
	Batch code		Single patient multiple use (for accessories only)
<b>IP21</b>	Protection against solid foreign objects and against harmful effects due to the ingress of water		

## WEEE- Waste Electrical and Electronic Equipment

### **THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.**

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

### **WARRANTY**

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions. The warranty shall not therefore apply in the case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national laws in the country of purchase, where appropriate.

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>		
The GCE858 is intended for use in the electromagnetic environment specified below. The customer or the user of the GCE858 should assure that it is used in such an environment.		
<b>Emission test</b> Emissions RF CISPR 11	<b>Compliance</b> Group 1	<b>Electromagnetic environment – guidance</b> The GCE858 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The GCE858 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>			
The GCE858 is intended for use in the electromagnetic environment specified below. The customer or the user of the GCE858 should assure that it is used in such an environment.			
<b>Immunity test</b> Electrostatic discharges IEC 61000-4-2	<b>IEC 60601 Test level</b> ±8 kV contact ±15 kV air	<b>Compliance level</b> ±8 kV contact ±15 kV air	<b>Electromagnetic environment - guidance</b> Floors should be wood, concrete or ceramic tile, synthetic material, the relative humidity should be at least 30 %.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>			
The GCE858 is intended for use in the electromagnetic environment specified below. The customer or the user of the GCE858 should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrical fast transient/burst IEC 61000-4-4	Power supply lines :±2 kV input/output lines: ±1 kV		Mains power quality should be that of a typical commercial or hospital environment.
Overvoltage IEC 61000-4-5	line(s) to line(s): ±1 kV line(s) to earth: ±2 kV 100 kHz repetition frequency		Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0,0% 300 cycle	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0,0% 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) Magnetic field IEC 61000-4-8	3 A/m 50 Hz/60 Hz	3 A/m 50 Hz/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

### Guidance and manufacturer's declaration - electromagnetic emissions

The GCE858 is intended for use in the electromagnetic environment specified below. The customer or the user of the GCE858 should assure that it is used in such an environment.

<b>Immunity test</b> Conducted RF IEC61000-4-6	<b>IEC 60601 Test level</b> 150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	<b>Compliance level</b> 150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	<b>Electromagnetic environment - guidance</b> Portable and mobile RF communications equipment should be used no closer to any part of the Smart including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. <b>Recommended separation distances</b> $d = 0,35 \sqrt{P}$ $d = 1,2 \sqrt{P}$
Radiated RF IEC61000-4-3	10V/m 80% Am at 1kHz	10V/m 80% Am at 1kHz	80 MHz a 800MH $d=1,2 \sqrt{P}$ Da 800 MHz a 2.7 GHz $d=2,3 \sqrt{P}$ where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretical y with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Guidance and manufacturer's declaration - electromagnetic emissions

The GCE858 is intended for use in the electromagnetic environment specified below. The customer or the user of the GCE858 should assure that it is used in such an environment

	Test Frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Modulation <sup>c)</sup> (W)	Distance (m)	Immunity Test Level (V/m)				
	385	380 – 390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1,8	0,3	27				
	450	380 – 390	GMRS 460, FRS 460	FM <sup>c)</sup> ± 5 kHz	2	0,3	28				
	710	704 – 787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217 Hz	0,2	0,3	9				
	745										
	780										
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0,3	28				
	870										
	930	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation <sup>b)</sup> 217 Hz	2	0,3	28				
	1 720										
	1 845										
	1 970										
	2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217 Hz	2	0,3	28				
	5 240	5 100 – 5 800						WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217 Hz	0,2	0,3
5 500											
5 785											

Note: If necessary to achieve the Immunity Test Level, the distance between the transmitting antenna and the ME Equipment or ME System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

<sup>a)</sup> For some services, only the uplink frequencies are included.

<sup>b)</sup> The carrier shall be modulated using a 50% duty cycle square wave signal.

<sup>c)</sup> As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The Manufacturer should consider reducing the minimum separation distance, based on Risk Management, and using higher Immunity Test Levels that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher Immunity Test Levels shall be calculated using the following equation:

$$E = 6/d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the Immunity Test Level in V/m

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GLOBALCARE MEDICAL TECHNOLOGY CO., LTD
Artsana S.p.A