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EN READ THESE INSTRUCTIONS AND KEEP FOR FUTURE REFERENCE

For your child's safety. WARNING!

⚠ WARNING! This is not a toy. Use under adult supervision. Throw away at the first sign of damage or weakness. 1.5V AA battery included.

• Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

• Cleaning and user maintenance shall not be made by children without supervision.

• Oral irrigators can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge under supervision.

• Children shall not play with the appliance.

REPLACING THE BATTERY

The battery can be replaced by unscrewing the bottom of the brush with the help of a tool. After replacing the battery, close the compartment and re-tightened the screw into position using the same tool.

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its

life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your retailer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the unwanted appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

The crossed bin symbol on the batteries or the product packaging indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your retailer when you buy similar rechargeable or non-rechargeable batteries.

Chemical symbols (Hg for mercury, Cd for cadmium, Pb for lead) printed beneath the crossed-out

wheeled bin symbol indicate the type of substance contained in the battery. The user is responsible for taking the batteries to special waste disposal sites at the end of their life so they can be processed and recycled. If dead batteries are collected correctly as separate waste they can be recycled, processed and disposed of ecologically; this avoids a negative impact on both the environment and human health, and contributes towards the recycling of the substances contained in the batteries. Improper disposal of the product by the user may harm the environment and human health.

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Artsana S.p.A. - Via Saldarini Catelli, 1